


Plastic Patterns

Plastic Preform Patterns

The highest quality materials are used to assure consistency and function. 8 cards per box

Each card — 2.25




Renfert Wax Patterns

Pre-Formed Wax Clasps

Green color, 10 per package.

Regular, Box of 15 — 9.99

Adhesive, Box of 15 — 9.99



Visit our online website at www.dentallabservices.com for specials, new products and used equipment listings!

Pattern Adhesive/Coats

Pattern Adhesive (Keystone)

A special liquid formulated to provide better adhesive characteristics for plastic preform patterns.

4 oz bottle — 8.45



Pattern Seal Spray (California Dental)

13 oz can — 9.20

Smooth & Stippled Casting Wax

Stippled Wax Sheets (Keystone)

Pressure sensitive. Green. Specify 22, 24 or 26 gauge.

Pkg of 15 — 25.25



Casting Wax (Renfert)

For use with chrome cobalt. Very easily adapted, tear-resistant, adheres well to the investment model. No elastic recovery. Specify .35, .40, .45, .50mm. Green.

Pkg of 15 — 9.99

Partial & Chrome Investment

Yellow-Stripe (Nobilium)

3-liquid, ethyl-silicate investment for chrome-cobalt partial dentures. Insures a smooth surface and dense grain structure.

2-50 lb. Powders + 1 Gallon of C+D with 470cc E Liquids — 110.00

50 lb. powder only — 24.00

1 gal. C+D w/470cc E Liquid — 67.00

X-20 Chrome (Whip Mix)

A phosphate investment for casting partials with chromium dental alloys. Mixes easily and produces a cast of unusual smoothness and strength. May be used with any good duplicating material intended for chrome investments.

P+L, 50 lbs — 68.00



Chrome Abrasives

Cut Off Wheels (Pacific Abrasives)

Non-contaminating, Italian made, aluminum oxide, ast cutting, reddish brown color.

(22*0.6mm) Box of 100 — 9.50

(25*0.7mm) Box of 100 — 12.25

(35*0.7mm) Box of 100 — 19.95

(40*1mm) Box of 100 — 23.00



Roughing Discs (Pacific Abrasives)

Non-contaminating, aluminum oxide, fast cutting.

(35*1.7mm) Box of 100 — 26.50

(32*2.2mm) Box of 100 — 26.50

(22*3.5mm) Box of 100 — 26.50

Chrome Welding Rods

Welding Rod, Chrome (CDP)

Easy flowing, high strength. Platinum color. Melting range 1700°-1800°F. Use with torch or electric.

10 dwt — 19.95

Duplicating Material

Caloid (California Dental)

1 gal — 27.50

3½ gal Pail — 75.00

Dupli-Coe-Loid (Coe)

Reversible duplicating hydrocolloid material for refractory models and pour denture resins.

1 gal — 57.00

Nobiloid "All Purpose" (Nobilium)

All purpose, reversible duplicating material for use with all investment types & can be used for duplicating stone models.

1 gal — 30.00

3½ gal — 89.00

Instaloid "Concentrated" (Nobilium)

All purpose, consistent, high accuracy for finest detail. Very versatile (small labs can mix "small" amounts). Just add water... instantly you have whatever volume of duplicating material you need.

1 gal (10 lb) — 77.00

95 oz (7 lb) — 57.00

Silicone Duplicating Material

Redu-It (American Dental)

Premium model duplicating silicone. Great dimensional stability and accuracy. Won't shrink on standing. Unlimited duplication can be made from master model at any time. Easy release from any model material.

2 Pints (A & B) — 74.00

2 Quarts (A & B) — 115.00

Electro Polishers

Electro-Polisher (California Dental)

Compact, solid state 10 amp unit. Polishes two cases at a time.

Unit — 425.00

Model Dip Hardener

Cal-Dip (California Dental)

Refractory model hardener. Use cold. Non-flammable. Safe for all investments.

1 gal Hardener — 39.50

Duplicating Flasks

Duplicating Flask (Ray Foster)

Flat Base — 18.00

Raised Base — 18.00



Pour Flask (DLS)

Each — 45.00

Chrome Polish

Cal-Strip (California Dental)

Chrome electrolytic solution, excellent for stripping high-fusing chrome alloys. For best results use at 125°F.

1 gal — 23.75

Diamonds

KIMBERLY DIAMONDS • ROBOT FG DIAMONDS • SINTERED • DIAMOND DISCS

Kimberly Diamond Instruments / Handpiece and Friction Grip


Kimberly Diamond Burs (Atwood)

Made in the USA. Proven quality since 1951. Specify Hand Piece or Friction Grip.

Pkg of 3 — 30.00

* FRICTION GRIP ONLY
† HAND PIECE ONLY

Pointed Tapers



Flat-End Tapers



Round-End Tapers



Rounds



Barrels



Footballs



Cones




Cylinders

Short Flat Heads



Long Flat Heads



Round Heads



Inverted Cones

Short Flat Heads



Long Flat Heads



Round Heads




Many Diamond Instruments are only available thru special order.

Shofu Diamond Robot Points

Friction Grip.

Pkg of 6 — 28.00



Sintered Diamonds

Sintered Diamonds (Keystone)

Hand piece grip burs.

Super Coarse Black: Chrome, Cobalt, Non-Precious and Hard Metals

Each — 33.95



Sintered Diamond Kit (Keystone)

Kit contains 12 various diamond burs. Burs are super coarse (black) hp-style.



Kit — 321.20

Diamond Cleaning Stone

Cleaning stone for sintered diamonds.

Each — 5.50

Flexis Diamond Discs (Keystone)



Flexible and versatile. For forming crown edges, For reducing crown borders to correct size before porcelain firing, For constructing and shaping of acrylic and ceramic surfaces, For infinitely variable comparison of boundary positions from metal to blending. For perfect finishing of the interdental area with bridges, Not to be used as a separating disc. Each box contains 8 discs plus mandrel.

Black (Coarse) — 41.00

Red (Medium) — 41.00

Green (Fine) — 41.00

White (X-Fine) — 41.00

Assortment — 41.00

Keystone Diamond Discs

Ti-Coated Diamond Discs

Titanium Nitrate Coated, sold each with one mandrel. Double sided, 22mm diameter.

SUPER PLASTER - For cutting plaster, ceramics and acrylics.



CONTOUR PERFORATED COARSE or FINE - For contouring ceramics and acrylics.

CUT FLEX COARSE - For cutting ceramics & acrylics and can be used for heavier cutting.



CUT FLEX FINE - For separating ceramics and acrylic

SUPER FLEX COARSE - Thin, flexible, for light contouring and working in interproximal slots, not cutting.

SUPER FLEX FINE - Single or Double sided, for separating ceramics, metals and acrylic.

SEP FLEX COARSE - For separating & contouring of ceramics, metals, plaster & acrylic.




TOP FLEX PERFORATED FINE - For contouring ceramics.

Single Pack — 16.50

Double Pack — 31.50

Horico Diamond Discs

Double sided: Diamond coating on both sides
Safe inside: Diamond coating on outside surface
Safe outside: Diamond coating on inside only
 All Diamond Discs listed below are 7/8" circular. Other sizes available or may be stocked in warehouse.



Diaflex (A)
 Full Coated 7/8" disc.

345-220 Double sided (.3mm thick) — **22.44**
 346-190 Safe Outside (.2mm thick) — **13.86**
 347-220 Safe Inside (.2mm thick) — **15.51**

Super Diaflex (B)
 Edge Coated 7/8" disc. Very-Flexible, Extra Fine.

355F-220 Double sided (.15mm) — **18.50**
 356F-220 Safe Outside (.1mm) — **17.50**
 357F-220 Safe Inside (.1mm) — **17.50**

Super Diaflex C (C)
 Edge Coated 7/8" disc. Very-Flexible, Med. Fine Grit.

355C-220 Double sided (.17mm) — **18.50**
 356C-220 Safe Outside (.12mm) — **17.50**
 357C-220 Safe Inside (.12mm) — **17.50**

Diaflex F (D)
 Full Coated 7/8" disc. Thin, Flexible, Extra Fine.

H358-220 Double sided (.18mm) — **22.50**
 H359-220 Safe Outside (.13mm) — **17.50**
 H360-220 Safe Inside (.13mm) — **17.50**

Super Diaflex S (E)
 Slotted 7/8" disc.

H368-220 Double sided (.18mm) — **22.50**

Perforated Diamond Discs

Double sided, 7/8" disc.

H350-220 (.3mm) — **22.50**

...Diagonal Slotted Super Diaflex

H353F-220 Flexible (.15mm) — **24.50**

H353C-220 Flexible (.17mm) — **24.50**

Model Prep Disc (Horico)

Double sided face, for model preparation.
 450mm (1-3/4") size.

Fine grit. (.3mm thick)

H365C-450 — **67.50**



DLS Diamond Discs

Diamond Disc

Double sided face, thin.
 (mandrel not incl.)

Each — **8.00**



George Taub Die Spacer

Tru-Fit Gold & Silver Die Spacers

A controlled technique for relieving dies before wax-up. Gives 25 microns of relief per coat. Use on stone, amalgam, epoxy and plated dies. Specify **Gold** or **Silver**.

½ oz bottle — **8.50**
 1 oz bottle — **15.00**
 2 oz bottle — **26.00**
 4 oz bottle — **42.00**



Space-It Red, Blue and Yellow

For those who prefer these colors rather than Tru-Fit gold or silver. These colors set up as fast as Tru-Fit, in less than 15 seconds. Each coat is 12-13 microns.

½ oz bottle — **9.50**
 1 oz bottle — **17.00**
 2 oz bottle — **28.00**
 4 oz bottle — **42.00**



Tru-Fit Thinner

Thinner for use with Tru-Fit and Space-It Die Spacers.

1 oz bottle — **4.50**
 2 oz bottle — **6.00**
 8 oz can — **13.00**
 16 oz can — **18.00**
 32 oz can — **25.00**

KerLab

Euro Classic Die Spacer

Manufactured to exacting tolerances in Germany, KerLab's reformulated line of die spacers, sealers and hardeners combine advanced chemistry with precision microtechnology to provide a consistently uniform coating with every application. Specify **Blue** or **Gray**.

½ oz bottle — **15.00**
 3.38 oz bottle — **45.00**



Euro Die Spacer Thinner

½ oz bottle — **15.00**
 100ml bottle — **45.00**

Euro Quick Set Die Spacer

Specify **Blue** or **Gray**.

½ oz bottle — **15.00**
 3.38 oz bottle — **45.00**

Euro Die Spacer Pens

Quick Set Die Preps available in a convenient Pen applicator with a polyfiber disposable tip.

3ml Blue SuperPen — **16.95**
 3ml Gray SuperPen — **16.95**
 3ml Die Hardener SuperPen — **16.95**



Tak

DieHard Die Spacer

Available in Red, Yellow, Blue, Orange or Clear.

4 oz bottle — **45.00**
 Thinner ½ oz — **7.25**



California Dental Products

Die Spacer

Clear, Blue, Orange, Silver or Gold.

½ oz bottle — **7.25**
 4 oz bottle — **37.00**
 Thinner 4 oz — **11.50**



Dental Ventures of America

Die Spacers

Widely used, smooth die spacer. Bonds chemically with DVA conditioner. Will NOT separate from the conditioner during waxing and/or finishing procedures. Available in Blue, Red, Silver or Gold.



1 oz Kit w/ Thinner — **19.00**

Die Lubricants

Die Lube (Neytech)

An excellent separating medium for wax patterns. It is non-oily and will not evaporate from the die, so it can be used sparingly. Since it is water soluble, it does not have to be washed off the pattern before investing. Acts as a surface tension reducer for the wax pattern.

½ oz bottle — **10.60**
 8 oz bottle — **16.40**
 16 oz bottle — **26.50**



Microfilm (KerLab)

A non-oily water soluble die lubricant that allows an accurate wax pattern to be removed from the die stone without distortion.

8 oz bottle — **18.50**
 32 oz bottle — **36.95**



Slikdie Lubricant (Slaycris)

Good for all types of crown & bridge-work. Especially useful in the reversible hydro-colloid technique. Water soluble lubricant can be rinsed off with cool tap water.

Slikdie 4 oz bottle — **12.10**



Die Lube (California Dental)

An excellent lubricant for stone, porcelain and acrylic. Wax patterns will easily release from any surface to which Die Lube has been applied.

4 oz bottle — **6.75**
 16 oz bottle — **20.25**
 32 oz bottle — **25.95**



Die Materials

DIE HARDENERS • CERAMIC SEPARATORS • DIE BLOCK-OUT • ADHESIVES

Die Hardeners

Die Model Hardener (KerrLab)

Best with KerrLab Die Spacers.

- Euro Classic Hardener, ½ oz — **15.00**
- Euro Classic Hardener, 3.3 oz — **45.00**
- Die Sealer/Hardener, ½ oz — **15.00**
- Die Sealer/Hardener, 3.3 oz — **45.00**

Stone Die Hardener (Taub)

Paints on in thin even coats, seals prosity, doesn't build up. Very cost effective.

- Resin, 2 oz — **12.00**
- Resin, 8 oz — **35.00**
- Resin, 16 oz — **55.00**
- Resin, 32 oz — **75.00**
- Thinner, 2 oz — **4.50**
- Thinner, 8 oz — **12.00**
- Thinner, 16 oz — **17.00**
- Thinner, 32 oz — **25.00**



Ceramic Separators

Cerama Seal & Sep (KerrLab)

Non contamination two liquid system allows building of porcelain directly to model without porcelain discoloration or sticking.

- Cerama Seal, ½ oz — **10.75**
- Cerama Sep, ½ oz — **12.50**



Rubber Sep (Taub)

A protective rubber die spacer, separator and release agent for temporary acrylic crowns and composite crown restorations. Specify Pink or White.

- ¾ oz — **15.00**
- 4 oz — **36.00**



Very Special Separator (DVA)

A die conditioner, which seals and hardens dies without adding space. Also, an oil-free die separator, which provides the most efficient separation of wax, porcelain or composite.

- Separator Pen, Pkg of 3 — **22.75**
- Separator, 8 oz — **38.00**
- Conditioner Pen, Pkg of 3 — **22.75**
- Conditioner, 8 oz — **38.00**
- Conditioner Thinner, 2 oz — **8.00**



Burn-Out Posts

Duralay Plastic Pins (Reliance)

For post and core. Serrated so Duralay will readily adhere. Tapered for easy insertion into the canal. Complete burnout with no residue.

- Plastic Pins (50) — **12.50**



Endowels (Star)

Endowel posts are tapered. Fast, most post preparations can be accomplished in less than 4 minutes. May be fitted with equal ease in round or elliptical canals. Available in Black, White, Yellow, Blue, Green or Assorted.

- Pkg of 50 — **27.00**



Die Marking Pencils

V.H.T. Marking Pencil (Whip Mix)

Mark margins on refractory dies. The blue color changes to red when heated and retains legibility through subsequent firing.

- Each — **7.50**



Die Block-Out

Zapit (Dental Ventures)

Plasticized Base Material with either spray or pen accelerator hardener.

- Base 1/3 oz — **13.75**
- Base 20gm — **29.50**
- Accelerator 2 oz — **8.95**
- Accelerator Pen (4) — **17.50**
- Accelerator 16 oz — **35.50**
- Solvent, 4 oz — **4.99**
- Tips for Base, Pkg of 10 — **10.50**



Adds It (ADDS)

- Base ½ oz — **11.50**
- Base 1 oz — **18.50**
- Accelerator 2 oz — **6.95**



Debubblers

Wax Pattern Cleaner (Jelenko)

A wetting agent for the immediate lowering of surface tension on wax patterns. Removes oil and decreases surface tension prior to investing.

- 8 oz spray — **12.00**



DeLar Surfactant (Delar)

Decreases contact angle of stone or plaster mix when poured in hydrocolloid, alginate, rubber silicone, wax or compound impressions.

- 8 oz spray — **19.95**
- 32 oz bottle — **47.95**



Bubble Free (California Dental)

A superior wax painting solution for smoother casting. Investment will flow uniformly and smoothly to all portions of the wax pattern.

- 16 oz bottle — **15.75**
- 32 oz bottle — **23.00**

Debubbler (KerrLab)

A water-base surface tension reducer that allows close adaptation of investment to all surfaces of the wax pattern. Apply by brushing, dipping or spray.

- 8 oz bottle — **16.50**
- 32 oz bottle — **32.50**



Smoothex (Whip Mix)

A surface tension reducing liquid for treating patterns to permit better adherence and flow of investment. Also excellent for treating synthetic impressions before pouring cast. The 4 and 8 oz bottles come with a fine spray pump.

- 8 oz bottle — **25.00**
- 16 oz bottle — **43.50**



Blue Debubbler (Mizzy)

Wax painting solution to produce accurate, smooth, bubble-free castings.

- 32 oz bottle — **14.40**
- Gallon Refill — **39.35**



Silicone Emulsion Spray Debubbler (Mizzy)

Produce smooth accurate, bubble-free surfaces with all impression materials and on all types of impressions.

- 4 oz bottle — **12.70**
- 16 oz bottle (Concentrated) — **38.70**

Adhesives

Insta-Cure Cynoacrylate (DLS)

- Super Thin, ½ oz — **4.50**
- Super Thin, 1 oz — **6.00**
- Gap Filling, ½ oz — **4.00**
- Gap Filling, 1 oz — **6.00**



Insta-Set Accelerator (DLS)

- 2.5 oz Spray Bottle — **4.50**
- 8 oz Refill — **12.00**

Zap CA (Pacer)

- Thin, Fast Setting ½ oz — **4.50**
- Thin, Fast Setting 1 oz — **7.50**



Zap-A-Gap (Pacer)

- Gap Filling ½ oz — **4.50**
- Gap Filling 1 oz — **7.50**

Zip Kicker Accelerator (Pacer)

- 2 oz Spray Bottle — **4.75**
- 8 oz Refill — **7.50**



Z-Ends (Pacer)

- Tips for Bottle
- Pkg of 10 — **4.50**

910 Adhesive (Permabond)

Thin fast setting adhesive.

- 3gm — **3.50**
- 1 oz — **11.50**



MDS QX-4 (ADDS)

Thin fast setting adhesive.

- 1 oz — **11.00**



Starbond (CPH)

Medium-set, high performance grade cyanocrylate.

- 1 oz — **8.50**



Accelerator Spray (Vivid)

Sets up fast.

- 420ml Spray Can — **15.00**