

Abrasives

SHOFU & DLS MOUNTED STONES




Green Mounted Stones (For Porcelain, Gold & Other Alloys)

Silicone Carbide

Dura-Green Mounted Stones (Shofu)

Made of high quality silicone carbide grits. Fast contouring and finishing of porcelain, composites, gold, silver and amalgam.

Pkg of 12 — 16.50




Green Mounted Stones (DLS)

Made from the finest quality silicon carbide with exclusive ceramic bond for the ultimate in rapid cutting and accuracy with minimal heat. Excellent for shaping, adjusting, finishing and fine grinding of Porcelain, Gold and other alloys.


Pkg of 10 — 4.95

Pkg of 100 — 19.95



Dura-White Mounted Stones (For Enamel, Composites & Porcelain)

Aluminum Oxide



Dura-White Mounted Stones (Shofu)

For fine finishing enamel, composites, and porcelain. Made of micrograined aluminum oxide grits, each stone is perfectly balanced to provide vibration-free performance.

Pkg of 12 — 15.25

White Mounted Stones (For semi & non-precious ceramic metal)

Aluminum Oxide

Lab-Series White Stones (Shofu)


For finishing semi & non-precious ceramic metal.

Pkg of 12 — 16.50

White Mounted Stones (DLS)

For semi and non-precious metals, for working without metal contamination.

Pkg of 10 — 4.95 Pkg of 100 — 19.95



Pink Mounted Stones (For precious ceramic metal)

Aluminum Oxide

Lab-Series Pink Stones (Shofu)


For finishing precious ceramic metal.

Pkg of 12 — 14.50

Pink Mounted Stones (DLS)

For nickel-chrome, cobalt, precious and semi-precious alloys, aluminum oxide, fine grain, non-contaminating.

Pkg of 10 — 4.95 Pkg of 100 — 19.95




Blue Mounted Stones (For partial denture and alum. oxide alloys)

Aluminum Oxide

Blue Mounted Stones (DLS)

For partial denture chrome-cobalt and aluminum oxide non-precious alloys.

Pkg of 10 — 4.95 Pkg of 100 — 19.95



Dedeco Grinding Wheels

7/8" x 1/8"

Fine Cut - For heatless chip-free finishing of porcelain, gold and enamel.

#5350 Box of 25 — **25.95**

Green Universals - Same qualities as fine cut, except "medium".


#5330 Box of 25 — **25.95**

Black Universals - For adjusting gold and metal castings.

#5320 Box of 25 — **25.95**

Soft Cut - Same qualities as fine cut, except "coarse".

#5310 Box of 25 — **25.95**



Heatless Wheels

Busch Silent Wheels (Pflingst)

Silicon carbide, for work on porcelain, fin and medium finish.

Box of 12 — **19.60**



13mm 16mm 19mm 22mm 25mm

2mm	761S	766S	771S	776S	781S
3mm	761S	766S	771S	776S	781S
4,5mm	761S	766S	771S	776S	781S

Heatless Wheels (Mizzy)

Grey Regular... Silicon carbide for the rapid reduction of plastic and metal materials.

White... Aluminum oxide & non-contaminating, for porcelain or metal.

Red... Aluminum oxide, non-contaminating, tough ceramic coating, use on all non-precious & semi-precious metals, porcelains, cercons & Lava & long lasting.

Regular, Box of 50 — **23.00**

White, Box of 50 — **23.00**

Red (sizes 1,3,5 or 8), Box of 50 — **23.00**



Wheel-	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
Dia-	1"	1"	1"	7/8"	7/8"	7/8"	3/4"	3/4"	3/4"	5/8"	5/8"	5/8"	1/2"	1/2"	1/2"
Width-	3/16"	1/8"	3/32"	3/16"	1/8"	3/32"	3/16"	1/8"	3/32"	3/16"	1/8"	3/32"	3/16"	1/8"	3/32"

"Koolies" Grinding Wheels (Dedeco)

Green... fine for porcelain and enamel.

Black... coarse - multi purpose favorite.

White... Non contaminating for finishing NP and semi-precious.

Green, Box of 50 — **34.95**

Black, Box of 50 — **34.95**

White, Box of 50 — **34.95**




Freedom Smooth Finishing Burs

Typhoon Sanding Sleeves

Leaves a smoother finish than other types of tungsten carbide burs. Typhoon burs last longer and do not load as quickly. Use the fine grits for general shaping and sanding. The coarse grit is for rapid and aggressive shaping and roughing out work.

*Call for pricing.



Sphere Taper Cylinder Ball-Nose Flame

Lathe Wheels

Lathe Wheels (Dedeco)

Green - Soft for grinding of porcelain

Black - Medium for grinding of porcelain

Green Medium - Medium for porcelain and metal.

Acrylic Base Plate Wheel - Coarse for base plates and acrylics

White Knock-Down Wheel - Coarsest for plaster, stone models and acrylic

3" x 1/4" Each — **8.95**

3" x 3/8" Each — **9.95**



"Jewelers" Lathe Wheel (Dedeco)

Soft and flexible for fine polishing of gold. Silicon carbide.

4" x 1/2" x 1/2" Each — **19.95**


Lathe Wheels (Wells)

For grinding porcelain teeth.

Choose from Green, Black or White.

3" x 1/4" Each — **7.25**

3" x 3/8" Each — **7.75**



Cratex Rubberized Lathe Wheels (Dixon)

Resists clogging, can be used wet or dry and cut freely smoothly.

Choose from Coarse Green, Medium Brown or Fine Red.

3" x 1/4" Each — **10.50**

3" x 3/8" Each — **12.50**

Green Grinding Wheels (Keystone)


For grinding porcelain teeth. Will not chip porcelain. 120 grit

3" x 1/4" Each — **7.90**

3" x 3/8" Each — **7.90**

2 1/2" x 3/8" Each — **6.40**


2 1/2" x 1/4" Each — **6.40**



Black Utility Grinding Wheels (Keystone)

For trimming denture base plate materials, stone models and acrylic.

3" x 1/4" Each — **7.40**



White Acrylic Knock-Down Wheels (Keystone)

For trimming base plates, acrylic and stone models. Coarse grit will not clog.

4" x 1/2" Each — **8.50**

3" x 3/8" Each — **7.25**



Abrasives

CUT-OFF WHEELS & DISCS • SEPARATING, PORCELAIN & PAPER DISCS

Pacific Abrasives Finishing Discs

Cut Off Wheels

Non-contaminating, Italian made, aluminum oxide, fast cutting, reddish brown color.

- (22•0.6mm) Box of 100 — **9.50**
- (25•0.7mm) Box of 100 — **12.25**
- (35•0.7mm) Box of 100 — **19.95**
- (40•1mm) Box of 100 — **23.00**



Roughing Discs

Non-contaminating, aluminum oxide, fast cutting.

- (35•1.7mm) Box of 100 — **26.50**
- (32•2.2mm) Box of 100 — **26.50**
- (22•3.5mm) Box of 100 — **26.50**

Dedeco Cutting Off Wheels

Cut-Off Wheels

Traditional favorites. Aluminum oxide for cutting and finishing chrome cobalt metals. Fast, cool, odor-free.

- (5163) 1¼" x .025", Pkg of 100 — **22.95**
- (5168) 1¼" x .032", Pkg of 100 — **22.95**
- (5166) 1¼" x .049", Pkg of 100 — **22.95**
- (5165) 1¼" x .062", Pkg of 100 — **22.95**
- (5169) 1½" x .032", Pkg of 100 — **22.95**
- (5167) 1½" x .062", Pkg of 100 — **22.95**

Roughing Discs

- (5164) 1¼" x .082", Pkg of 100 — **32.95**

Rubi "Slitters"

For non-precious and semi-precious metals.

- (5174) (38•0.8mm), Pkg of 100 — **32.95**

"Red" Discs

Faster cutting and finishing of chrome cobalt metals.

- (7022) 1¼" x .062", Pkg of 100 — **22.95**
- (7023) 1½" x .025", Pkg of 100 — **22.95**
- (7024) 1½" x .040", Pkg of 100 — **22.95**

"Fastcut" Discs

For both chrome cobalt and non-precious metals.

- (5211) (25•0.7mm) Pkg of 100 — **15.95**
- (5212) (35•0.7mm) Pkg of 100 — **27.95**
- (5213) (35•0.6mm) Pkg of 100 — **32.95**
- (5214) (40•1mm) Pkg of 100 — **27.95**

California Dental Products

Nylo-Disks

Non-contaminating for alloys. Aluminum Oxide.

- Regular 1" x 0.7" Box of 100 — **21.75**
- Thin 7/8" x .06" Box of 100 — **21.75**

Nylo-Cutters

Non-contaminating for alloys. Aluminum Oxide.

- 1-1/4" x 3/32" Box of 25 — **18.25**

Keystone Separating Discs

Joe Dandy Discs

Universal use for fast cool cutting on porcelain, gold, acrylic, alloys and others.

Specify 1/2", 5/8", 3/4", 7/8" or 1" disc. .023" thick.

- Flat Double-Cutting, Box of 100 — **10.70**
- Flat Safe Side, Box of 100 — **10.70**
- Cup Double-Cutting, Box of 100 — **20.70**



Cut-Off Wheels

Non-contaminating, aluminum oxide, reddish brown.

- 1¼" x 0.62" Box of 100 — **21.50**
- 1½" x 0.40" Box of 100 — **21.50**
- 1½" x 0.25" Box of 100 — **21.50**

Roughing Discs

- (35•1.7mm) Box of 100 — **26.50**
- (32•2.2mm) Box of 100 — **26.50**

Red Flash Discs

For all chrome, precious and semi-precious alloys. Aluminum oxide - 25,000rpm.

- 15/16" x .025" Box of 100 — **10.90**



Veri-Thin Discs

Suited for cutting and finishing porcelain and gold.

Silicon carbide - Grey - 10,000-20,000rpm.

- 7/8" x 0.15" Box of 100 — **18.50**
- 3/4" x 0.15" Box of 100 — **18.50**

Dura-Thin Discs

Non-contaminating companion to the ultra-thin disc. Made of aluminum oxide and formulated for fine cutting and finishing of ceramic metals.

- 7/8" x .010" Box of 100 — **13.70**

Dedeco Separating Discs

"SILICONE CARBIDE"

Wafer Thin Discs

Very thin discs for porcelain, gold, enamel, acrylic and crown & bridge metals. Pack of 100.

- (5178) Double side, 7/8" x .015" — **34.95**
- (5180) Safe side, 7/8" x .015" — **34.95**

Stripper Discs

"Coarse" black discs for fastest finishing and cutting of porcelain and gold. Pack of 100.

- (5177) 7/8" x .025" — **19.95**

"ALUMINUM OXIDE" A/O Thin Discs

Ultra-thin end cutting discs - for gold and all metals.

- (5173) 7/8" x .009" Box of 25 — **11.95**

Rubi Mini Slitters

For fastest, cool cutting of all non-precious & semi-precious ceramic metals. Pack of 100.

- (5176) Mini Slitters, 7/8" — **17.95**

High Speed Discs

For fastest, cool cutting of gold and all metals. 100 pcs.

- (5209) (22•0.6mm) — **9.95**

Renfert Separating Discs

For Precious Metals

- #71 Ultra-Fine, X-thin 22•0.17 (50) — **46.50**
- #72 Thin, non-flexible 22•0.3 (100) — **27.50**
- #74 Universal 22•0.3 (100) — **32.50**

For Non-Precious Metals

- #57 Standard 25•06 (100) — **11.50**
- #68 Standard 32•1.6 (100) — **35.00**
- #66 Standard 37•1.0 (100) — **32.50**
- #65 Standard 38•0.6 (100) — **32.50**

For Porcelain

- #75 Thin 22•0.38 (50) — **22.50**
- #60 Extra-Thin 22•0.2 (50) — **18.00**
- #74 Universal 22•0.3 (100) — **32.50**

Porcelain Discs

Ultra-Thin Disc (Keystone)

The industry standard for fine cutting and finishing on porcelain interproximals and margins.

Grey 7/8" x .009"

- Box of 25 — **13.70**



Speed-Eze (California Dental)

For porcelain, gold and chrome alloys.

- 7/8" Double Sided Box of 100 — **21.50**
- 7/8" Single Sided Box of 100 — **21.50**

NM Slims (Dedeco)

Ultra-thin end cutting discs. For porcelain, gold, enamel and crown & bridge metals. Non-contaminating.

- 7/8" x .009" Box of 25 — **11.95**

Ultra-.009's (Dedeco)

Ultra-thin end cutting discs. For porcelain, gold, enamel and crown & bridge metals.

- 7/8" x .009" Box of 25 — **11.95**

Speedy Discs (Union Broach)

The original porcelain disk.

- 7/8" Flat Double-Sided Box of 100 — **22.50**
- 7/8" Flat Single-Sided Box of 100 — **22.50**
- 7/8" Cup Double-Sided Box of 100 — **22.50**

Moore Paper Discs

Paper Discs, Brass Center

Choose from 1/2", 5/8", 3/4" or 7/8" diameter.

- EMERY X-Fine, Fine, Med, Coarse, X-Coarse
- GARNET X-Fine, Fine, Med, Coarse, X-Coarse
- SAND X-Fine, Fine, Med, Coarse, X-Coarse
- Small Pkg/50 — **2.95**



MOORE

Dedeco Mounted Rubber Polishers

Midgets

Color coded and long-lasting mounted rubberized abrasives for pinpoint anatomical finishing and polishing of all materials.



- BROWN:** Long lasting for pre-finishing all metals.
 - BLACK:** Coarsest for amalgam and gold.
 - GREEN:** Traditional for fast polishing of amalgam and gold.
 - BLUE:** Hi-shine for final polishing of amalgam, precious & semi-precious metals.
 - IMPS(RED):** Coarse for removal of tartar and tobacco stains. Eliminates need of pumice.
 - WHITE:** Ultra-fine for highest polish on composites, enamel and gold.
 - GREY:** For instant hi-glaze touch-up polishing of porcelain.
- HP Shape, Box of 12 — **20.95**

Mounted Silicone Polishers

Brownie, Greenie, Supergreenie Points (Shofu)

Made of silicone impregnated with high quality polishing ingredients. For fast brilliant finishes on amalgam, gold and cast alloys.

- Brownie:** achieve smooth surfaces
- Greenie:** for a lustrous polish
- Supergreenie:** for a superpolish



HP or FG Box of 12 — **21.95**

Acrylic Polishers

Jumbos (Dedeco)

Rubberized abrasives are designed for heavy duty polishing and finishing of acrylics and metals. HP shank.



- BLACK:** Coarse (silicon carbide) for finishing metal, trim acrylics.
 - GREEN:** Medium (silicon carbide) for polishing metals, smooth acrylics.
 - WHITE:** Very fine (alum. oxide) for highest polish on metals and acrylics.
- Box of 5 — **18.95**

Arbor Bands & Abrasive Rolls

Arbor Bands (Keystone)

Sizes available are 1/4", 3/8", 1/2" or 3/4". Choose from Coarse, Medium or Fine.



- Bands, Bag of 100 — **13.50**
- Bands, Box of 500 — **86.65**

Heavy Bonded Arbor Bands (Keystone)

Long lasting, coarse grit. Sizes available: 3/8" or 3/4".


- Heavy Bonded Bands, Bag of 100 — **13.25**
- Heavy Bonded Bands, Bag of 500 — **53.75**

Keystone Silicone Polishers

Deluxe Silicone Wheels & Points

BLACK MEDIUM - For all metals.
WHITE COARSE - For all metals, amalgam, porcelain and composites.
BLUE FINE - For final polish of porcelain or metals.
PINK COARSE - For porcelain.
 Packaged 100 per box.

- Black Square Edge — **55.00**
- Black Knife Edge — **55.00**
- Black Points — **49.50**
- White Square Edge — **55.00**
- White Knife Edge — **55.00**
- White Points — **49.50**
- Blue Square Edge — **55.00**
- Blue Knife Edge — **55.00**
- Blue Points — **49.50**
- Pink Square Edge — **55.00**



Keystone Superior Rubber Wheels and Points

DARK GREY REGULAR - For excellent finish on chrome or polish on gold.
DARK GREY HARD - For an excellent finish on chrome and polish on gold.
LIGHT GREY MED-SOFT - For fast scratch-free finishing of chrome.
BLUE GLOSS - For a lustrous finish on chrome, semi-precious and precious metals.
BLUE FLEX - Softer than Blue for both chrome and noble metals.
GREEN MED-FINE - All purpose for all metals. Fast high quality finishing.
BLACK MED-COARSE - Hard and flexible for fast chrome finishing.
RED FINE - Provides a high quality finish on precious and non-precious.
WHITE - Non-contaminating for ceramic gold and porcelain - Aluminum Oxide.
GREEN UTILITY COARSE - For fast, effective polishing on gold, porcelain and gold.

		Sq. Edge	Knife Edge	Points	
SUPERIOR RUBBER POLISHING WHEELS AND CLASP POLISHERS		Dk Grey R	12.85	----	13.10
		Dk Grey H	12.85	----	13.10
		Lt Grey	12.85	----	----
		Blue Gloss	12.85	23.40	13.10
		Blue Flex	12.85	----	13.10
		Green MF	12.85	23.40	13.10
		Black MC	12.85	23.40	13.10
		Red F	12.85	23.40	13.10
		White	17.45	----	----
		Green UTIL	15.75	23.40	----
7/8" x 1/8"					
Box of 100					

Dedeco Rubberized Favorites

BROWN - For finishing & polishing of chrome cobalt, non & semi-precious metals.
BLACK - Hard for finishing and polishing chrome cobalt metals.
GREEN COARSE - For scratch removal & polishing of chrome cobalt metals.
BLUE - High glaze for polishing non-precious, semi-precious & chrome cobalt.
RED FLEXIES FINE - Soft and flexible for fine polishing of gold.
WHITE FLEXIES - Soft and flexible for ultra fine polishing of porcelain & gold. NC.
GREY - Medium hard for finishing & polishing chrome cobalt metals.
GREEN MEDIUM - For fast finishing of all metals.

		Sq. Edge	Knife Edge	Points	
SUPERIOR RUBBER POLISHING WHEELS AND CLASP POLISHERS		Brown	26.95	30.95	26.95
		Black	18.95	22.95	18.95
		Green C	18.95	22.95	18.95
		Blue	18.95	22.95	18.95
		Red Flex	22.95	26.95	22.95
		White Flex	24.95	28.95	----
		Grey	18.95	----	18.95
		Green M	22.95	----	----
	7/8" x 1/8"				
	Box of 100				

Burlew and Cratex Rubber Wheels

Burlew Rubber Wheels (Jelenko)

Pleasantly peppermint flavored for use in the mouth.

- Wheel 5/8" Box of 50 — **16.00**
- Sulci 3/8" Box of 50 — **12.00**

Cratex Rubber Wheels

Choose: Green-Coarse, Brown-Medium, Red-Fine or Light Green-Xfine

- Wheels Box of 100 — **29.75**
- Points Box of 50 — **17.50**
- Cylinders Box of 50 — **17.50**

Keystone Truing Stone

Point & Wheel Truer

A silicon carbide block for dressing and cleaning mounted abrasive points. Reshaping points back to perfect concentricity is done by running the mounted point in the groove comparable to size.



Each — **6.15**

Abrasives

SILICONE WHEELS & POINTS • HATHO DISCS • MANDRELS • ECONO CUTTERS • BLASTING MATERIAL

Pacific Silicone Wheels/Points

Pre-Polish Polishers

For precious & semi-precious alloys.



Bag of 10 — **4.50**

Box of 100 — **42.60**

Dark Grey - Soft for excellent scratch-free finish.

G-1 Square Edge 7/8"

G-2 Knife Edge 7/8"

G-3 Cylinder 7/8" x 1/4"

White - Flexible for excellent scratch-free finish.

G-A Square Edge 7/8"

G-B Knife Edge 7/8"

Porcelain Touch-Up Polishers

For precious and semi-precious alloys. For touch-up on porcelain, laminate veneers.

Bag of 10 — **4.95**

Box of 100 — **50.75**

SG-1 Coarse (Grey) Sq. Edge 7/8"

SG-1K Coarse (Grey) Knf. Edge 7/8"

SG-1C Coarse (Grey) Cylinder 7/8"

SG-3 Medium (Lt. Blue) Sq. Edge 7/8"

SG-3K Medium (Lt. Blue) Knf. Edge 7/8"

SG-3C Medium (Lt. Blue) Cylinder 7/8"

SG-4 High-Shine (Lt Grey) Sq. Edge 7/8"

SG-4K High-Shine (Lt Grey) Knf. Edge 7/8"

SG-4C High-Shine (Lt Grey) Cylinder 7/8"

SG-1M Imp. Coarse (Grey) Sq. Edge 7/8"

SG-1K1M Imp. Coarse (Grey) Knf. Eg. 7/8"

Hi-Shine, Gold Flash Wheels

Pink color. For precious and semi-precious alloys.

For beautiful high shine without any rouge.

Square Edge (GF4) or Knife Edge (GF5).

Box of 100 — **54.00**

N.P. Special

For non-precious alloys. For very smooth finish.

Square Edge or Knife Edge available.

Box of 100 — **42.60**

Mini Points

Choose from Coarse, Medium or Hi-Shine.

Box of 100 — **48.60**

Mandrel — **4.30**

Hatho Discs

Habras Discs (Hatho)

A three-dimensional abrasive/flexible design that provides a uniform finish. Ideal for heavy cleaning (coarse), removing light oxide (medium), pre-polishing, semi-finishing, blending, texturing and polishing (fine). Removes oxides and polish without using compounds. Ideal for cleaning metal work and hard-to-reach areas. Available in Yellow-80grit, Red-220grit, Blue-400grit, Pink-Pumice, Orange and Green.

Pkg of 4 with Mandrel — **6.70**

Pkg of 4 — **5.35**

Pkg of 48 — **60.25**



Habras Disc Kit — **49.40**

Dedeco Diamond Sanding Caps

"Dia-Dome" for the Lathe

Lathe sanding caps designed to replace arbor bands and traditional acrylic trimmers. Works with standard 3/4" arbor band mandrels on 1/4" shanks.

(2741) Coarse, Pkg of 2 — **54.95**

(2742) Fine, Pkg of 2 — **54.95**



"Dia-Dome" for the Lathe

Uniquely designed for shaping and contouring.

(2745) 3/32" shank, Pkg of 2 — **3.95**

(2744) 3/32" shank, 1/2" dia., Pkg of 2 — **16.95**

Mandrels / HP shank

Arbor Band Mandrels (Keystone)

Choose from 1/4", 3/8", 1/2" or 3/4" head.

Mandrel — **2.45**

Rubber Head — **.70**



Handpiece Mandrel

For separating discs. 3/32"

Pkg. of 12 — **9.50**



Tapered HP Mandrel

For separating discs, tapered style. 3/32"

Pkg. of 12 — **9.50**



Large Head HP Mandrel

For thin discs, 8mm head. 3/32"

Pkg. of 12 — **12.50**



Spiral HP Mandrel

For rubber points and bullets. 3/32"

Spiral Mandrel Pkg. of 12 — **9.50**

Spiral Mandrel w/ nut Pkg. of 12 — **9.50**



Corkscrew HP Mandrel

For felt points and felt cylinders. 3/32"

Pkg. of 12 — **9.50**



Moore's Mandrels

Snap-On Mandrel

For Moore's brass center discs.

#7 HP Each — **4.00**



Econo-Cutters

"Fibre-Cut" Discs (Dedeco)

Shatterproof, 1/32" thick, reinforced discs cut everything from hard metals to soft acrylics.

(7001) 1 1/2" Each — **1.75**

(7002) 2 1/2" Each — **2.75**

(7003) 3" Each — **3.25**

(7004) 4" Each — **4.25**



"Fibre-Cut" Mandrel (Dedeco)

Gives extra durability and finer performance for discs 7002, 7003 & 7004.

Each — **4.95**




Econo-Cutters (Keystone)

Fast cutting and long lasting. Use at speeds up to 24,000 rpm.

1 1/2" Each — **2.75**

2 1/2" Each — **3.50**

4" Each — **4.95**



Econo-Cutter

Mandrel (Keystone)

Special mandrel needed for 2-1/2" and 4" Econo-Cutter discs.

Each — **5.40**



Blasting Compounds

Aluminum Oxide (DLS)

White 25 Micron 13 lb — **29.50**

White 25 Micron 52 lb — **89.50**

White 50 Micron 13 lb — **26.00**

White 50 Micron 52 lb — **85.00**

Brown 50 Micron 13 lb — **26.00**

Brown 50 Micron 52 lb — **85.00**

White 100 Micron 13 lb — **26.00**

White 100 Micron 52 lb — **85.00**



Glass Beads (DLS)

50 Micron 13 lb — **18.50**

50 Micron 52 lb — **60.50**

Walnut Shellblast (California Dental)

Abrasive for shallow indentations. Suitable for the removal of oxides and for the roughening porcelain surfaces before firing.

20 lb — **43.80**

Cobra Blast (Renfert)

Abrasive for shallow indentations. Suitable for the removal of oxides and for roughening porcelain surfaces before firing.

20 lb — **43.80**

