



King Hue cross linked teeth


Bioform shades, New Hue Molds.

Each Card... 1.75

King Hue Shade Guide... 5.00

*Buy 100+ cards @ only 1.50 per card

Shades: 50,59,62,65,66,67,69,77,81,87*



KingHue Anteriors - Upper & Lower

	A	B	C	
B3	8.30 - 6.70 - 38.00			
	8.30 - 4.50 - 31.00			
3M	9.00 - 7.00 - 38.50			
	8.20 - 5.00 - 31.50			
2D	9.50 - 7.8 - 42.00			
	9.20 - 5.00 - 32.00			
3D	9.30 - 8.00 - 44.00			
	8.80 - 5.00 - 34.00			
A24	9.50 - 8.00 - 43.50			
	8.80 - 5.00 - 34.00			
A25	9.50 - 8.50 - 46.00			
	9.50 - 5.00 - 34.00			
A26	9.80 - 9.00 - 47.00			
	10.00 - 6.00 - 38.00			
2N	10.30 - 7.50 - 42.00			
	9.00 - 5.00 - 33.00			
3N	10.00 - 8.00 - 45.00			
	9.50 - 5.00 - 34.00			

2P	11.30 - 8.50 - 47.00	
	9.80 - 5.50 - 36.00	
3P	10.50 - 8.50 - 48.00	
	10.00 - 6.00 - 38.00	
1H	12.00 - 8.50 - 48.00	
	10.80 - 5.80 - 38.50	
4H	10.30 - 8.80 - 48.00	
	9.30 - 5.50 - 37.50	
133	8.80 - 7.50 - 41.00	
	8.30 - 4.50 - 31.00	
135	10.00 - 8.80 - 47.50	
	10.00 - 6.00 - 38.00	
263	10.00 - 7.50 - 41.00	
	9.00 - 5.00 - 33.00	
264	10.80 - 8.00 - 45.00	
	10.00 - 6.00 - 38.00	
266	12.00 - 9.00 - 49.00	
	11.50 - 5.50 - 38.00	
267	12.40 - 9.50 - 51.00	
	11.50 - 5.50 - 38.00	

KingHue Posteriors - Upper & Lower

31Z 0°	
X 29.50 - U 7.20	
Y 31.00 - L 8.40	
32F 10°	
X 31.00 - U 8.50	
Y 32.50 - L 9.00	
33F 10°	
X 32.00 - U 9.00	
Y 34.50 - L 9.40	
29M 20°	
X 29.00 - U 7.20	
Y 32.00 - L 7.80	
29L 20°	
X 29.00 - U 8.20	
Y 32.50 - L 8.20	
31M 20°	
X 31.00 - U 7.90	
Y 32.50 - L 8.50	
31L 20°	
X 30.90 - U 8.70	
Y 34.50 - L 8.70	
33M 20°	
X 32.00 - U 8.60	
Y 35.50 - L 8.70	
30M 33°	
X 29.00 - U 7.70	
Y 31.50 - L 7.70	
30L 33°	
X 29.00 - U 8.00	
Y 31.00 - L 8.40	
32M 33°	
X 31.00 - U 7.60	
Y 33.50 - L 8.20	
32L 33°	
X 30.00 - U 8.40	
Y 33.00 - L 8.70	
34M 33°	
X 32.50 - U 8.80	
Y 34.50 - L 9.30	
34L 33°	
X 33.00 - U 8.90	
Y 35.50 - L 9.20	

Articulation	POSTERIOURS			
	0°	10°	20°	33°
B3	-	-	29M	30M
3M	-	-	29M	30M
2D	-	32F	29L	30L
3D	31Z	33F	31M	30M
A24	-	-	29M	30M
A25	31Z	33F	31M	30M
A26	31Z	33F	31M	30M
2N	-	-	29L	30L
3N	31Z	32F	31M	30L
2P	31Z	33F	31L	30L
3P	-	33F	33M	34M
1H	-	33F	-	34L
4H	-	33F	33M	34M
133	-	-	29M	30M
135	31Z	33F	33M	32M
263	-	32F	29L	30M
264	31Z	-	31L	30M
266	-	-	33M	34L
267	-	-	33L	34L

Posteriors	
0°	10° 20° 33°

Teeth / Ultrasonics

VERATONE MULTI CHROMATIC TEETH • L&R ULTRASONICS • SONIX ULTRASONICS




Veratone multi chromatic cross linked acrylic teeth

Bioform shades, New Hue Molds.

Each Card... **3.75**

Veratone Shade Guide... **10.00**

Shades: **A1,A2,A3,A4,B1,B2,B3,C2,D3**




Veratone Anteriors - Upper & Lower

13	9.1 - 7.8 - 41.3 8.0 - 5.1 - 33.1	
25	9.6 - 8.0 - 41.7 8.2 - 5.2 - 33.0	
35	9.7 - 8.0 - 43.1 8.6 - 5.5 - 35.5	
36	9.9 - 7.9 - 42.9 8.4 - 5.4 - 35.3	
56	9.7 - 8.2 - 43.9 8.6 - 5.5 - 35.5	
57	10.8 - 8.4 - 45.4 9.3 - 5.1 - 34.4	
66	10.0 - 8.0 - 42.9 8.6 - 5.4 - 35.4	
76	10.2 - 8.1 - 44.8 9.3 - 5.7 - 35.9	
77	10.3 - 8.5 - 45.3 9.7 - 5.6 - 38.0	

78	10.0 - 8.3 - 45.5 9.0 - 5.7 - 36.1	
85	9.5 - 8.5 - 46.3 8.6 - 5.6 - 35.9	
88	10.5 - 8.6 - 46.1 8.6 - 5.6 - 36.2	
98	11.3 - 8.9 - 47.5 11.7 - 5.6 - 37.4	
99	12.5 - 9.7 - 51.4 10.6 - 6.1 - 41.4	

Articulation

Posteriors



anteriors	posteriors
13	41 12
25	41 11
35	41 44
36	41 11
56	12
57	44 13
66	12
76	41 11
77	44 13
78	41
85	44 13
88	44 13
98	45
99	45

KingHue Posteriors - Upper & Lower

11	A 29.5 B 31.3 C 8.1 D 7.1	
12	A 30.2 B 32.0 C 9.0 D 7.7	
13	A 30.4 B 33.4 C 8.5 D 7.35	
41	A 29.2 B 31.0 C 8.8 D 8.3	
44	A 32.0 B 34.0 C 9.1 D 8.6	
45	A 34.5 B 35.9 C 9.4 D 8.4	

L&R Ultrasonics

Quantrex units have black vinyl panel and cover. Sweepzone units are constructed with AgION antimicrobial coated stainless steel wrappers.

Model PC3 - 18 oz. Tank

White, self contained, compact, high-performance. Each — **155.00**

Quantrex Q90 - 2 Quart

Exterior: 6½ x 6½ x 8½ - Interior: 5½ x 5½ x 4
With Timer — **325.00**
With Timer & Heat — **375.00**

Quantrex Q140 3¾ Quart

Exterior: 10½ x 6½ x 8
Interior: 5½ x 5½ x 4
With Timer & Drain — **375.00**
With Timer, Heat & Drain — **450.00**



*****With Sweepzone Technology*****

Sweepzone Ag200, 4¼ Quart

With Timer & Drain — **525.00**
With Timer, Heat & Drain — **575.00**

Accessories

600ml Beaker — **14.50**
600ml Beaker Top — **9.75**
For Q140 and S140 units



L&R Ultrasonic Solutions

General Purpose Cleaner

Non-Ammoniated, Gallon — **29.00**
Ammoniated, Gallon — **29.00**

Temporary Cement Remover

Advanced Formula, Gallon — **33.00**

Plaster & Stone Remover

Gallon — **31.00**



Tartar, Lt. Stain, Perm. Cement Remover

Gallon — **29.00**

Sonix IV

Sonix SS-Series Cleaners

Powerful cavitation, automatically adjusts liquid levels, temperatures and solutions for varying work loads.

Sonix SS116 - 2 Quart

Exterior: 6½ x 5½ x 8½ - Interior: 5½ x 5 x 4
Each — **235.00**
With Heater & Basket — **260.00**

Sonix SS117 - 3 Quart

Exterior: 10 x 5½ x 9.1 - Interior: 9½ x 5½ x 4
Each — **320.00** With Heater — **335.00**

Sonix SS117 - 6.4 Quart

Exterior: 12 x 6.4 x 11.6 - Interior: 11½ x 5½ x 6
Each — **502.00**

Dentsply Neytech Ultrasonics

Mini ULTRAsonik

Mini ultrasonic cleaner, 8 oz. capacity, 5¼ x 5¼.

Each — 75.00



Ney 14H

1.5 quart, timer, cover and heater. 7.5" x 3.5" x 3"

Each — 250.00

Ney 19H

2 quart, timer, cover and heater. 6" x 5.5" x 4"

Each — 260.00

Ney 28H

3 quart, timer, cover and heater. 11.6" x 6" x 6"

Each — 360.00

Heraeus Kulzer Solutions

VK-1 Gen. Purpose Cleaner

Gallon — 36.00

VK-3 Plaster & Stone Remover

Gallon — 32.00

VK-4 Tartar & Stone Remover

Gallon — 32.00

VK-6 Temp. Cement Remover

Gallon — 32.00

Super Concentrated Solutions

Ultradose (L&R)

Rapidly cleans blood, debris, and other matter from instruments. Super concentrated: 1 oz. = 1 gallon solution.

24-1 oz Packets — 34.00



Ultradose Gen. Purpose Cleaner

Super concentrated: 33gm packet = 1 gallon solution.

24-33 gm Packets — 34.00



Ultradose Denture Cleaner Powder

Removes tartar and light stains in ultrasonic. Safe for cleaning acrylic, porcelain, and most metals. 5gm packet = makes 400ml

24-5 gm Packets — 28.00



Neyclean GP (Neytech) Cleaner & Degreaser

Makes 75 gallons of solution.

Each — 38.70

Neyclean EZ (Neytech) Enzymatic Cleaner & Presoak

Makes 75 gallons of solution.

Each — 44.50

Neyclean SP (Neytech) Stone, Plaster & Inv. Remover

Each — 22.30

Sonic Power Granules (CDP)

Ultrasonic Cleaning Granules. General purpose, non-ammoniated cleaning. Bio-degradable, with rust inhibitor. Suitable for removing dried blood, polishing compound, sand cleaning of prosthetics.

16 oz (Makes 32 gallons) — 25.50

Porcelain Stripping Solutions

EZ-Strip

Removes oxidation, investment and porcelain from precious and non-precious alloys.

16 oz Bottle — 8.00

Gallon — 26.50

Super Stripper

Industrial strength hydrofluoric acid substitute. Cleans castings and strips porcelain better than any other available.

Kit (1 pint Stripper + Neutralizer) — 18.50

2 Pint — 27.50

4 Pint — 43.50

Silicone Sprays

All purpose release agent for use on metals, rubber, waxes, plastics, epoxies, plaster, stone, amalgam and compounds.

California Dental Products

16 oz Can — 10.50

Lang Dental

16 oz Can — 12.50

Dental Resources

16 oz Can — 15.00



Vacuum Forming Machines



The Machine III

New easy latch system, strongest 2-stage motor, large sample pack of materials, all aluminum construction and a 5-Year Warranty.

T&S — 285.00



Air Vac XQ

Finest vacuum former available, European style processing, easily attaches to your air line for the finest appliance and bleaching trays, 5-Year Warranty.

Mizzy — 287.00

Accu-Vac (Buffalo)

Uses compressed air-powered suction, creates highly accurate vacuum adaptations, gasket-free frame accepts 5" diameter round & 5" x 5" square materials.

Each — 275.00



Sta-Vac II (Buffalo)

Versatile new conversion table combination, vacuum table and bead basket, high performance electric vacuum forming system. With 2 ConverTable tables.

Each — 525.00



Econo-Vac (Buffalo)

Illuminated rocker switches alert operator if heater is on, lifetime silicone gasket never needs replacing, square balanced base will not tip, 2-Year Warranty.

Each — 325.00



Proform 1 (DR)

Hexagonal post that keeps sliding frame from rotating. All aluminum construction to eliminate rust. No gasket in plastic holder. 5-Year Warranty.

Each — 285.00



Proform 2 (DR)

Dual chambered, heavy duty motor for greater suction output. Spring loaded frame to accommodate a variety of materials. 5-Year Warranty.

Each — 350.00



Pressure Dome (DR)

For fabricating better fitting conventional dental appliances when using vacuum forming techniques. This attachment hooks into any existing compressor air line (needs 80 PSI line pressure with 3/8" line) to provide extra pressure when adapting thermoplastics over an existing model. It will fit any existing type 5" x 5" vacuum machine.

Each — 95.00



Proform E Gasket

Proform E Gasket (DR)

Process a gasket denture start-to-finish within 24 hours.

Clear, Pkg of 12 — 47.50

Pink, Pkg of 12 — 47.50

Vacuum Forming / Wax

VACUUM FORMING MATERIAL • BASEPLATE, BITE, BOX, CARVING & DIAGNOSTIC WAXES

Pro-Form Vacuum Forming Materials (Keystone Industries)



Pro-form Soft E.V.A.

Soft EVA Clear Mouthguard

A soft versatile material used in bleaching, model duplication, soft trays, etc. Soft, clear, EVA.

- .040 EVA, Pkg of 25 — 15.00
Also great for Bleaching Trays
- .060 EVA, Pkg of 25 — 17.50
- .080 EVA, Pkg of 25 — 17.75
- .120 EVA, Pkg of 25 — 20.00
- .120 Blue, Pkg of 25 — 20.00
- .160 Clear, Pkg of 25 — 20.00
- .160 Blue, Pkg of 25 — 20.00
- .200 Clear, Pkg of 25 — 23.00
- .200 Blue, Pkg of 25 — 23.00

Pro-form Splint/Tray Resins

Temporary Splint/Surgical Tray

An excellent material for splints, surgical trays, retainers, and partials. Clear color.

- .020, Pkg of 25 — 5.50
- .020, Pkg of 50 — 11.00
- .040, Pkg of 25 — 5.75
- .040, Pkg of 50 — 17.00
- .060, Pkg of 25 — 15.00
- .060, Pkg of 50 — 27.50
- .080, Pkg of 25 — 16.00
- .080, Pkg of 50 — 30.00
- .100, Pkg of 25 — 18.75
- .100, Pkg of 50 — 36.00

Pro-form Base Plate & Tray Resins

Base Plate & Tray Material

Used to make base plates and custom trays. Custom trays include spacer foam.

- .060 Base Plate Pink, Pkg of 25 — 7.00
- .080 Base Plate Pink, Pkg of 25 — 8.50
- .100 Base Plate Pink, Pkg of 25 — 9.75
- .125 Custom Tray Blue, Pkg of 25 — 11.50
- .150 Custom Tray Blue, Pkg of 25 — 13.50

Pro-form Specialty Sheet Resins

Crown & Bridge (Temporaries)

Used as a matrix for making temporaries.

- .020 Clear, Pkg of 25 — 11.50

Coping Material

Great for temporaries, soft niteguards, and retainers.

- .020 Clear, Pkg of 25 — 6.50
- .030 Clear, Pkg of 25 — 7.50
- .040 Clear, Pkg of 25 — 11.00

Retainer Material

A hard, thin material used in the production of invisible retainers. Bonds to acrylic.

- .030 Clear, Pkg of 25 — 16.50

Base Plate Waxes

Premium Base Plate Wax (Corning)

- Pink, 5 lb — 32.50

Modern Pink #3 (Heraeus Kulzer)

Soft, flows well.

- Pink, 1 lb — 18.50
- Pink, 5 lb — 65.00

Base Plate Wax

(Keystone)

General purpose, economical wax.

- Pink, 5 lb — 52.00
- Red, 5 lb — 52.00



Nobilwax (Nobilium)

5 lb. box

- #1 Pink, Extra Hard — 29.00
- #2 Pink, All Season-Thick — 29.00
- #2 Pink, All Season-Thin — 34.50

Base Plate Wax (Yates & Bird)

Pink is excellent wax, simulates flesh tones.

- All Season Wax, Pink 5 lb — 38.50
- Red Refined Wax, 5 lb — 38.50

Bite Wax

Bite Wafers (Corning)

U-shaped wax gives quick, precise check bites and occlusal transfers.

- Box of 56 — 13.80
- Box of 192 — 38.70

Bite Blocks

No waste, ready to use.


- Pink 3/8", 112 per box — 16.60
- Bite Sticks, 100 per box — 16.60

Bite Blocks (Modern Materials)

- Pink soft, 120/pkg — 54.50
- Pink hard, 120/pkg — 54.50

Bite Blocks (Keystone)

- Pink 3/8" — 27.50



Bite Blocks (Keystone)

- Pink — 23.50

GEO Anatomic Wax (Renfert)

Anatomically shaped wax bite rims, for occlusal surfaces, time savings during modeling, wax-ups and temporary models can be created quickly.

- Box — 68.50



Riji Form Bite Rims (Yates & Bird)

A unique ridge-formed bite rim that saves time, saves wax, saves labor. No waste. Shape to the model, cut off and use exactly what you need.

- 4 lb Box — 27.00

Boxing/Strip Wax

Boxing Wax (Corning)

- Regular Strip Wax Red, 1 lb — 9.50
- Thin Strip Wax Red, 1 lb — 9.50
- Thin Strip Wax Red, 5 lb — 32.50

Modern Materials (Heraeus Kulzer)

- Regular Red, 1 lb — 15.25
- X-Thin Red, 1 lb — 19.25

Boxing Wax Strips and Sticks (KerrLab)


- Sticks 6 oz, 175g — 00.00
- Strips 1 lb, 454g — 00.00



Carving Wax

GEO Classic Modeling Wax (Renfert)

For work using the bunsen burner flame. Colors available: Beige Transp, Blue Transp, Mint Transp, Beige Opq, Gray Opq, Mint Opq, or Tooth color opaque.



- 75g — 21.88

(Whip Mix) Life-Like Presentation Wax

- 12 sticks, bx — 20.25
- Covered tin, 42g — 20.25



Casting Wax Sheets

Casting Wax Sheets (Keystone)

- | PINK | GREEN |
|------------------|------------------|
| 20 gauge — 15.00 | 20 gauge — 15.00 |
| 22 gauge — 15.00 | 22 gauge — 15.00 |
| 24 gauge — 15.00 | 24 gauge — 15.00 |
| 26 gauge — 15.00 | 26 gauge — 15.00 |
| 28 gauge — 15.00 | 28 gauge — 15.00 |
| 30 gauge — 15.00 | 30 gauge — 15.00 |

Plastic Stipple Sheets (Nobilium)

- 22 gauge, medium grit — 4.50
 - 22 gauge, coarse grit — 4.50
 - 24 gauge, medium grit — 4.50
 - 24 gauge, coarse grit — 4.50
 - 26 gauge, medium grit — 4.50
 - 26 gauge, coarse grit — 4.50
- GREEN COLOR

Diagnostic Wax

Diagnostic Wax (Caliber)

- Snow White, 75 gm — 10.50
- Tan, 75 gm — 10.50

Diagnostic Wax (D-R Miner)

Available in Snow White or Ivory.

- Box of 20 sticks — 26.50
- 6 oz Tin — 26.50

Dipping Wax

Magic Dipping Wax (Bego)

Changes color when ready to carve.

4 oz Bottle — **32.50**



"Art-Ese" Granular Dipping Wax (Caliber)

Violet, Orange, Yellow or Green.

6 oz Bottle — **15.50**



Corning Dip Wax (Corning)

Comes in bead form, for wax copings.

Green 8 oz — **24.00**

Yellow 8 oz — **24.00**

Dip 'N Cast Wax Dipping Wax

(Dental Ventures of America)

The unique composition does not require that a technician "reseat" the margins in most cases. The reduction of elimination of the "resealing" procedure can cut 60% to 75% of the production time of a coping.

Colors: Flexilbe Green, Flexible Orange, Firm Green or Firm Yellow.

1/2 lb — **53.00**

1/4 lb — **32.00**

Dipping Wax (KerrLab)

Cubit form, available in Green, Yellow, Blue or Red.

9 oz Bottle — **17.95**

18 oz Bottle — **24.95**



GEO Dipping Wax (Renfert)

Systematic wax copings for processing within an electronic wax pot. 200g

Yellow — **24.00**

Orange — **24.00**

Green — **24.00**



Pico Dip Wax for hotty LED

Yellow — **24.00**

Orange — **24.00**

Gator Dipping Wax (Whip Mix)

Minimal distortion/shrinkage. 88 wafers/box - 8 oz.

Dip Wax Yellow — **48.50**

Build-Up Wax Blue — **48.50**



Glue Gun

Hotstikz (Dental Ventures)

Specially sticky wax used with glue gun applicator for easy cleanup.

Blue Sticks, 5 lb — **83.75**

Glue Gun only — **39.50**



Inlay-Casting Wax

Inlay Wax Regular (Corning)

Green Sticks, 1 lb — **31.30**

Regular, Hard or Special Hard.

Blue Sticks, 1 lb — **31.30**

Regular, Hard or Special Hard.

Green Chunks, — **18.50**

Regular, Hard or Special Hard.

Blue Chunks, 1 lb — **18.50**

Regular, Hard or Special Hard.

Maves Inlay Wax

Scientifically prepared wax with correct and reliable working qualities.

Green Sticks, 1 lb

#1 Medium Soft — **38.00**

#2 Stiff — **38.00**

#3 Hard Pliable — **38.00**



Green Chunks, 1 lb

#1 Medium Soft — **28.00**

#2 Stiff — **28.00**

#3 Hard Pliable — **28.00**

Ivory Chunks, 1 lb

#2 Stiff — **28.00**

Inlay Casting Wax (KerrLab)

Blue Regular

12 sticks — **10.50**

Blue Regular

120 sticks — **32.95**

Blue Hard, 12 sticks — **11.50**

Blue Hard, 120 sticks — **32.95**

Green Hard, 1 lb chunk — **24.50**

Green Hard, 120 sticks — **35.50**

Ivory Regular, 12 sticks — **11.50**



Bellewax Sculpture Wax

(KerrLab)

2 oz Tin — **12.50**

Green, Blue, Red, Yellow, Brown, Green X-Hard, Bio

Pontic Blue, Dentin, Translucent or Pontic Green.

9 oz Bottle (wax cubits) — **23.50**

Green, Blue, Red, Yellow, Brown, Green X-Hard, BioPontic Blue or Dentin.

IQ 4⁺ Sculpturing Wax

(Keystone/Yeti)

Higher inner stability, very hard, not fragile, higher point of solidification of 64° C - the ideal wax for summer, unique rigidity, very good sculpturing features, easy to carve and does not contain synthetic compounds.

45 gm Tin — **21.50**

Blue, Beige, Neon Green, Grey, Neon Yellow, Ivory Diagnostic or Black Cervical.



Thowax Sculpturing Wax

Highest quality European inlay wax.

60 gm Tin — **21.00**

Grey, Beige, Blue, Maroon, Violet, Green, Beige-Maroon, Cervical Red, Bordeaux or Milling Green.

Slaycris Wax (Slaycris Labs)

Regular, 3 oz Tin — **8.50**

Med-Soft wax especially for use in conjunction with occlusal buildup, formulated for a reversible hydrocolloid bridge and inlay technique.

Hard, 3 oz Tin — **8.50**

To be used with wax #1 for warmer climates.



Creative Color Sets

Chromowax Set (Benzer)

Shade Creation Set — **125.00**

Contains 10 differently colored waxes, 1 felt pad.

Esthetic Wax

(KerrLab)

For a naturally colorful appearance in diagnostic wax-ups and training in ceramic colorization.

Esthetic Wax Kit — **33.75**

Contains 2 oz. Ivory color tin, 2 oz. Tin of assorted wax colors, 1 belle wax duster brush.



Margin Wax

Margin Wax (Bego)

Blue Medium, 2 oz — **32.50**

OTM Occlusal Waxes (Corning)

OPAQUE GRAY, HARD - Has an optical color, allowing the anatomical carving to be seen with ease.

Tin — **13.75**

TRANSLUCENT DARK BLUE, HARD - Has an optical color, allowing the anatomical carving to be seen with ease.

Tin — **13.75**

TRANSLUCENT DARK GREEN, HARD - A carving wax created to sculpt accurate thickness in a coping or framework design. The translucency allows the tech to determine degree of depth.

Tin — **13.75**



OTM Carving Waxes

OPAQUE BEIGE, HARD - A versatile wax, can be used as an occlusal or carving wax. the opaque quality prevents eye strain.

Tin — **13.75**

OPAQUE LIGHT BLUE, HARD - Excellent wax for sculpturing, its high strength allows for easy carving without chipping, flaking or crumbling.

Tin — **13.75**

OPAQUE LIGHT GREEN, HARD - This is softer for carveable knife-edge margin finishing.

Tin — **13.75**

OTM Margin Wax

DARK RED - Unlike the others, this wax has an excellent ability to retain its integrity to create a perfect margin that flows easily and adheres to other waxes.

Tin — **13.75**



Wax

BLOCKOUT, UNDERCUT, PRESSABLE, ADHESIVE, STICKY, UTILITY & WIRE WAXES

Thowax Underlining Wax

Smooth flowing, a new dimension of precision.

70 gm Tin, Bordeaux — **26.50**



Margin Wax (KerrLab)

Burgundy, 2 oz Tin — **12.50**

Yellow, 2 oz Tin — **12.50**

Margin Wax (MDL)

Blue, 2 oz Tin — **11.00**

Blue, ½ lb — **23.50**

Red Margin & Base Wax (Slaycris)

Produces easily visible margin edges and accurately draws off over imperfections in the die to precisely define the anatomy in metal restorations.

Red, 3 oz Tin — **12.50**

Blockout/Undercut Waxes

Block Out Compound (Buffalo)

Each — **11.00**

"Art-Ese"

Blockout Wax (Caliber)

Yellow 75 gm Tin — **10.50**



GEO Block-Out Wax (Renfert)

Optimum bond on plaster or metal when blocking out and preparing models for duplication. Its high opacity means it can be seen even in thin layers. Its high melting point of 86°C (185°F) prevents softening when duplicating. Pink/Opaque

75g — **15.04**

GEO Under-Cut Wax

Dead wax for small undercuts and cavity coverage before inlay modeling. Excellent for flooding when forming secondary parts made of wax. It is softer and more elastic than cervical wax.

75g — **11.82**

Undercut Wax (Neytech)

24 gm Tin — **9.50**



Blockout/Undercut Waxes

Carving Wax Blocks (D-R Miner)

Available in 5/8" or 1/2". Choose Ivory, Green, Pink, Blue or Red blocks.

Pkg of 20 — **25.00**

Starter Wax Kit

Starter Kit Wax (Hi-Tech)

Kit — **45.00** Contains six different types of wax.



Pressable Wax

Pressable Ceramic Wax (HiTech)

Opaque wax designed for full crown pressable ceramics. This wax is specially designed to have clean burnout without any residue and/or color contamination. Choose Beige, Turquoise, White or Emerald color.

90 gm Tin — **29.50**

Pressable Margin Wax (HiTech)

Yellow, 60 gm Tin — **25.75**

Pure White Wax (Slaycris)

An ash free wax created especially for complete residue-free burnout.

White, 3 oz Tin — **14.50**

Pressure Sensitive Wax

Adhesive Wax (Keystone)

Flexible 4" x 4" wax sheets, coated on one side. Available in 20, 22, 24, 26, 28 or 30 Gauge.

Box of 32 sheets — **20.50**

Sticky Wax

Sticky Wax (KerrLab)

Orange Sticks
Pkg of 12 — **10.95**

Orange Sticks
Pkg of 120 — **32.95**



Sticky Wax (Corning)

Yellow Lumps, 1 lb — **13.30**
Yellow Sticks, 1 lb — **28.30**
Red Lumps, 1 lb — **15.30**
Red Sticks, 1 lb — **28.50**

Thowax Sticky Wax (Yeti)

Beige, 60gm tin — **23.75**



Sticky Wax (MDL)

Available in Yellow, Pink or Blue (sticks only).

Chunks, ½ lb — **18.00**
Sticks, ½ lb — **23.00**
Sticks, 1 lb — **43.00**



Sticky Wax (D-R Miner)

Orange Sticks — **22.00**
Natural Sticks — **22.00**

Utility Wax Strips

Utility Wax (Keystone)

For built up and post-damming purposes.

Red Round Strips, 1 lb — **9.50**
White Round Strips, 1 lb — **9.50**
Red Square Rope, 1 lb — **9.50**
White Square Rope, 1 lb — **9.50**



Utility Wax (Hygenic)

Square White Ropes, 1 lb (Box of 44) — **11.50**
Square Red Ropes, 1 lb (Box of 44) — **11.50**
Round White Strips, 1 lb (Box of 80) — **9.50**

Utility Wax (Corning)

Round Red, 1 lb — **8.50**

Wax Solvent

Wax Solvent #618 (CDP)

Safe, efficient, pleasant scent. For wax elimination and wax dressing of acrylic cases.

Gallon — **29.95**

Wire Wax

Ready Made Wire Shapes (Kerr)

Round - Choose 8,10,12,14,16,18 or 20 Gauge.

Pkg of 28 (4" long) — **15.95**



Wax Wire Spools (Kerr)

1/2 lb - Choose 6, 8,10,12,14,16,18 or 20 Gauge.

½ lb Roundspool — **24.50**

Kewax Wire Wax (Keystone)

Choose 6, 8,10,12,14,15,16,18,20,21,22 or 24 Gauge.

½ lb Spool — **20.00**

Kewax ½ Round Wire Wax

Choose 6,8,10,12 or 14 Gauge.

½ Round, Spool — **20.00**

Wire Wax Spools (Corning)

Choose 6, 8,10,12,14,16,18, or 20 Gauge.

½ lb Spool — **20.00**

Geo Wax Wire (Renfert)

Tension-free, won't revert back to its original shape, can be individually matched to object due to the varying diameters, easy to shape but remains dimensionally stable, no knotting of ripping during shaping.

Available in Blue Medium-Hard or Turquoise Hard.

2.0 mm, 250g — **12.00**

2.5 mm, 250g — **12.00**

3.0 mm, 250g — **12.00**

3.5 mm, 250g — **12.00**

4.0 mm, 250g — **12.00**

5.0 mm, 250g — **12.00**

